

# Acu•Dispo•Cautery® / Acu•Cautery®

Disposable and Reusable Surgical Cauteries and Replacement Tips

## Acu•Dispo•Cautery®

Disposable Surgical Cautery

Loop Tip 2200°F

Fine Tip 2200°F

2" Fine Tip 2200°F

Fine Tip 1300°F

Fine Tip 850°F

## Acu•Cautery®

Reusable Surgical Cautery

Fine Tip 2200°F

Loop Tip 2200°F

Fine Tip 1300°F

Made in the U.S.A.



**Acuderm® inc.**  
Real Service/Real Value



Intertek

CMDCAS








1-800-327-0015

# Surgical Cauteries from Acuderm



Acuderm's disposable and reusable cauteries are available in high and low temperature settings, plus tip styles designed for specific surgical applications.

## Disposable Cauteries

 Acu-Dispo-Cautery® FT 2200°F, Fine Tip STERILE. <b>CE</b>	Fingertip control, Safety cap. Uses 2 non-replaceable AA alkaline batteries. STERILE.	DCFT
 Acu-Dispo-Cautery® LT 2200°F, Loop Tip STERILE. <b>CE</b>	Fingertip control, Safety cap. Uses 2 non-replaceable AA alkaline batteries. STERILE.	DCLT
 Acu-Dispo-Cautery® FT2 2200°F, 2" Fine Tip STERILE. <b>CE</b>	Fingertip control, Safety cap. Uses 2 non-replaceable AA alkaline batteries. STERILE.	DCFT2
 Acu-Dispo-Cautery® M1 Small size cautery 1300°F, Fine Tip STERILE. <b>CE</b>	Fingertip control, Safety cap. Uses 1 non-replaceable AA alkaline battery. STERILE.	DCM1
 Acu-Dispo-Cautery® MM1 Mini size cautery 850°F, Fine Tip STERILE. <b>CE</b>	Fingertip control, Safety cap. Uses 1 non-replaceable AAA alkaline battery. STERILE.	DCMM1

## Cautery Kits

Acu-Cautery® Kit 2200°F or 1300° F Temp 2 Fine, 2 Loop Tips incl. <b>CE</b>	Reusable cautery with 4 replaceable tips. Incl. 2 replaceable AA alkaline batteries. NON-STERILE.	ACK (2200°F) ACK LT (1300°F)
---	---	---------------------------------



# Acuderm® inc.

5370 NW 35 Terrace, Suite 106

Fort Lauderdale, FL 33309

**Toll-Free: (800) 327-0015**

Tel: (954) 733-6935 • Fax: (954) 486-3602

E-mail: [cust-service@acuderm.com](mailto:cust-service@acuderm.com)

Web: [www.acuderm.com](http://www.acuderm.com)

Acuderm products meet applicable world-class standards, including

**Quality Certified and Registered ISO 13485:2003**

**Certified to Medical Device Directive CE**

**CMDCAS, UL 60601-1**




## For Reliable, Pinpoint Hemostasis

Choose Acuderm surgical cauteries for quick pinpoint hemostasis. Choose the convenience of our disposable Acu-Dispo-Cautery® or the cost-effectiveness of our reusable Acu-Cautery®, Acuderm gives you a choice of high (2200°F) or low (1300°F) temperature cautery.

Lightweight and easily maneuverable, these cauteries are ideal for use in a wide range of neurosurgery, ophthalmic and cosmetic surgery procedures. Replacement tips for Acu-Cautery® (our reusable cautery) are readily available in fine or loop tips, plus specialty tips such as an angled tip, vasectomy tip or gyn tip.






**Safe for use on patients with pacemakers.**

## Reusable Cauteries

 Acu-Cautery® 2200°F FT 2200°F, Fine Tip. <b>CE</b>	Fingertip control, Safety cap. Uses 2 replaceable AA alkaline batteries. Incl. 1 Fine Tip. NON-STERILE	ACFT
 Acu-Cautery® 2200°F LT 2200°F, Loop Tip. <b>CE</b>	Fingertip control, Safety cap. Uses 2 replaceable AA alkaline batteries. Incl. 1 Loop Tip. NON-STERILE	ACLT
 Acu-Cautery® 1300°F FT 1300°F, Fine Tip. <b>CE</b>	Fingertip control, Safety cap. Uses 1 replaceable AA alkaline battery. Incl. 1 Fine Tip. NON-STERILE	FT 1300

## Replacement Tips

NOTE: Replacement tips will fit any Acu-Cautery® (2200°F FT, 2200°F LT, 1300°F FT models). All replacement tips are STERILE, high temp (2200°F) and come with handle covers.

 Angular Fine Tip	 Fine Tip	Fine Tip	ACT-FT
 Vasectomy Tip	 Loop Tip	Loop Tip	ACT-LT
 5" GYN Loop Tip		Angular Fine Tip	ACT CA
		Vasectomy Tip	ACT VT
		5" GYN Loop Tip	ACT 5LT

Distributed By: