

Acu•Dispo•Curette®

Sterile Disposable Dermal Curette



Patented fit and feel for
precise curettage.
Made in U.S.A



Acuderm® inc.
Real Service/Real Value



Intertek


CMDCAS




1-800-327-0015

Acu•Dispo•Curette®

Available In Six Sizes:




A Stainless Steel, Disposable Curette Available In 2 Tip Styles:



Cup

1mm cup for precise placement and hard-to-reach areas such as nailbeds.




Loop

Loop tip is available in 5 sizes, 2mm to 7mm.

- Acu•Dispo•Curette® is consistently sharp every time.
- Acu•Dispo•Curette® is ergonomically designed for that “fit and feel” required for dermal curettage. No other disposable curette has this patented design.
- Try a sampler pack - five each of every size loop-end curette (stock no. RCOM) or five each of every size, including loop-end and cup currettes (stock no. RCOMPL).
- For “made-in-America” quality, look for the white handle and the Acuderm® logo.

Stock Numbers

Size	Box of 25	Box of 50	Box of 100
1mm	R0125	R0150	R0100
2mm	R0225	R0250	R0200
3mm	R0325	R0350	R0300
4mm	R0425	R0450	R0400
5mm	R0525	R0550	R0500
7mm	R0725	R0750	R0700
All Sizes		RCOMPL	Box of 30
5 Loop Sizes		RCOM	Box of 25



Each Acu•Dispo•Curette® is packaged in a sterile peel-apart pouch. These are available in boxes of 25, 50 or 100.

Distributed By:



5370 NW 35 Terrace, Suite 106

Fort Lauderdale, FL 33309

Toll-Free: (800) 327-0015

Tel: (954) 733-6935 • Fax: (954) 486-3602

E-mail: cust-service@acuderm.com

Web: www.acuderm.com

Acuderm products meet applicable world-class standards, including

Quality Certified and Registered ISO 13485:2003

Certified to Medical Device Directive CE

CMDCAS, UL 60601-1