

Acu•Punch® Kit

A Complete Disposable Kit for Skin Biopsy Procedures



Acu•Punch® Kit PLUS
Includes Webster Needleholder 4^{3/4}"

Acu•Punch® Kit ULTRA
Includes Webster Needleholder 4^{3/4}"
and Acu-Suture® 10" Nylon, C-2, 5/0

Made in U.S.A.



Acuderm® inc.
Real Service/Real Value



Intertek

CMDCAS



1-800-327-0015

Acu•Punch® Kit

Available In 3 Standard Choices or Customize.

A Complete Disposable Kit for Skin Biopsy Procedures

All the instruments are sterile and ready to use.

- Perfect for use at nursing homes or on hospital consults.
- Ideal for HIV and other high risk patients.
- Eliminate cross-contamination risks
- Cost-effective: eliminates sterilization and prep time.

Stock Numbers

| Description | Packaging | Cat # |
|--|-----------|------------|
| Acu-Punch Kit with 3mm Punch | Box of 20 | PK320 |
| Acu-Punch Kit with 4mm Punch | Box of 20 | PK420 |
| Acu-Punch Kit Plus with 3mm Punch | Box of 20 | PK320 PLUS |
| Acu-Punch Kit Plus with 4mm Punch | Box of 20 | PK420 PLUS |
| Acu-Punch Kit Ultra with 3mm Punch | Box of 20 | PKU320 |
| Acu-Punch Kit Ultra with 4mm Punch | Box of 20 | PKU420 |

Customize the Acu-Punch Kit with your choice of 14 punch sizes, sutures, and forceps (plastic or metal).



Acu•Punch® Kit

Each kit contains:

- 1 Acu•Punch® (3mm or 4mm)
- 1 Acu-Scissors (3-1/2" curved)
- 2 Gauze pads
- 1 Acu-Forceps (disposable)
- 1 Povidone iodine prep pad
- 1 Alcohol prep pad
- 1 Adhesive bandage (7/8" round)
- 1 CSR wrap (12" X 12")

Acu•Punch® Kit Plus

Each kit contains:

- 1 Acu•Punch® Kit (3mm or 4mm)
(See details above)
- Plus
- 1 Webster Needle Holder (4-3/4")



Acu•Punch® Kit Ultra

Each kit contains:

- 1 Acu•Punch® Kit (3mm or 4mm)
(See details above)
- Plus
- 1 Webster Needle Holder (4-3/4")
- 1 Acu-Suture® (10" black nylon, C-2 needle, 5-0 suture)



5370 NW 35 Terrace, Suite 106
Fort Lauderdale, FL 33309

Toll-Free: (800) 327-0015

Tel: (954) 733-6935 • Fax: (954) 486-3602

E-mail: cust-service@acuderm.com

Web: www.acuderm.com

Acuderm products meet applicable world-class standards, including

Quality Certified and Registered ISO 13485:2003

Certified to Medical Device Directive CE

CMDCAS, UL 60601-1

Distributed By: